



Kompatibilitätsübersicht der Fisherbrand SureOne™ Filterspitzen

Diese Liste erhebt keinen Anspruch auf Vollständigkeit. Sollten Sie Ihre Pipette in dieser Liste nicht finden, kontaktieren Sie uns bitte unter der Tel. 07851-99457-0 oder per e-mail : sales@novodirect.de

Pipetter Description	5984D	5984C	5984F	5984H	5984L	5984N
BIOHIT, , eLine, e10, ch: MC8, vol:0.2-10	yes	yes				
BIOHIT, , eLine, e10, ch: SC, E, vol:0.2-10	yes	yes				
BIOHIT, , eLine, e300, ch: MC8 E, vol:10-300						
BIOHIT, , mLINE®, m10, ch: MC12, vol:0.5-10	yes	yes				
BIOHIT, , mLINE®, m10, ch: MC8, vol:0.5-10	yes	yes				
BIOHIT, , mLINE®, m10, ch: SC, vol:0.5-10	yes	yes				
BIOHIT, , mLINE®, m100, ch: SC, vol:10-100					yes	
BIOHIT, , mLINE®, m1000, ch: SC, vol:100-1000						yes
BIOHIT, , mLINE®, m200, ch: SC, vol:20-200						
BIOHIT, , mLINE®, m300, ch: MC12, vol:30-300						
BIOHIT, , mLINE®, m300, ch: MC8, vol:30-300						
BIOHIT, , PROLINE® ELECTRONIC, , ch: MC-8, E, vol:25-250					yes	
BIOHIT, , PROLINE® ELECTRONIC, , ch: MC-8, E, vol:5-100					yes	
BIOHIT, , PROLINE® ELECTRONIC, , ch: SC, E, vol:5-100					yes	
BIOHIT, , PROLINE®, , ch: SC, vol:5-50					yes	
BIOHIT, , PROLINE®, PROLINE®, ch: SC, vol:100-1000						yes
BIOHIT, , PROLINE®, PROLINE®, ch: SC, vol:20-200					yes	



Pipetter Description	5984D	5984C	5984F	5984H	5984L	5984N
BIOZYM, , PRECISION, , ch: SC, E, vol:0.2-10	yes	yes				
BRANDTECH, , TRANSFERPETTE®, , ch: MC-12, vol:10-100						
BRANDTECH, , TRANSFERPETTE®, , ch: MC-12, vol:5-50						
BRANDTECH, , TRANSFERPETTE®, , ch: SC, vol:0.5-10	yes	yes	yes	yes		
BRANDTECH, , TRANSFERPETTE®, , ch: SC, vol:10-100						
BRANDTECH, , TRANSFERPETTE®, , ch: SC, vol:20-200					yes	
BRANDTECH, , TRANSFERPETTE®, , ch: SC, vol:5-50						
BRINKMANN, , TRANSFERPETTE®, , ch: MC-12, vol:5-25						
BRINKMANN, , TRANSFERPETTE®, , ch: MC-8, vol:50-200					yes	
CappAero, , , CappAero, ch: MC-8, vol:25-200						
CappAero, , , CappAero, ch: SC, vol:0.5-10	yes	yes				
CappAero, , , CappAero, ch: SC, vol:25-200					yes	
CappAero, , , Softline, ch: MC8, vol:50-300						
CappAero, , , Softline, ch: MC8, vol:5-50						
CappAero, , , Softline, ch: SC, vol:0.5-10	yes	yes				
CappAero, , , Softline, ch: SC, vol:10-100					yes	
CappAero, , , Softline, ch: SC, vol:25-200					yes	
CappAero, , , Softline, ch: SC, vol:5-50					yes	
CappAero, , , Softline, ch: SC, vol:5-50						
COSTAR®, , 8 PETTE, , ch: MC8, vol:20-200					yes	
COSTAR®, , ONEPETTE™, , ch: SC, vol:10-100						
COSTAR®, , ONEPETTE™, , ch: SC, vol:1-20						
COSTAR®, , ONEPETTE™, , ch: SC, vol:20-200					yes	
EPPENDORF®,REFERENCE, Physio Care, ch:SC,vol:100-1000						yes
EPPENDORF®,REFERENCE, Physio Care, ch: SC,vol:50-200					yes	



EPPENDORF®,REFERENCE, Physio Care/Metal Cone,ch:SC,vol:2-20	yes	yes	yes	yes		
Pipetter Description	5984D	5984C	5984F	5984H	5984L	5984N
EPPENDORF®,REFERENCE,Physio Care/Metal Cone, ch:SC, vol:5-10	yes	yes	yes	yes		
EPPENDORF®,RESEARCH, Physio Care, ch: MC 12, vol:10-100					yes	
EPPENDORF®,RESEARCH, Physio Care, ch: SC, vol:0.1-2.5	yes					
EPPENDORF®,RESEARCH, Physio Care, ch:SC,vol:100-1000						yes
EPPENDORF®,RESEARCH, Physio Care, ch:SC,vol:10-100					yes	
EPPENDORF®,RESEARCH,Physio Care, ch:SC, vol:20-200					yes	
EPPENDORF@, , RESEARCH, Physio Care, ch: SC, vol:2-20					yes	
EPPENDORF@, 2000, REFERENCE, , ch: SC, vol:0.1-2.5	yes					
EPPENDORF@, 2000, REFERENCE, , ch: SC, vol:0.5-10	yes	yes	yes	yes		
EPPENDORF@, 2000, REFERENCE, , ch: SC, vol:10-100					yes	
EPPENDORF@, 2000, REFERENCE, ch: SC, vol:200-1000						yes
EPPENDORF@, 2000, REFERENCE, , ch: SC, vol:50-200					yes	
EPPENDORF@, 2000, REFERENCE, GRAY, ch: SC, vol:2-20	yes	yes	yes	yes		
EPPENDORF@, 2000, REFERENCE, YELLOW, ch: SC, vol:2-20					yes	
EPPENDORF@, 2100, RESEARCH, , ch: MC-8, vol:0.5-10			yes			
EPPENDORF@, 2100, RESEARCH, , ch: MC-8, vol:10-100					yes	
EPPENDORF@, 2100, RESEARCH, , ch: SC, vol:0.1-2.5	yes					
EPPENDORF@, 2100, RESEARCH, , ch: SC, vol:100-1000						yes
EPPENDORF@, 2100, RESEARCH, , ch: SC, vol:10-100					yes	
EPPENDORF@, 2100, RESEARCH, , ch: SC, vol:20-200					yes	
EPPENDORF@, 2100, RESEARCH, , ch: SC, vol:2-20						
EPPENDORF@, 2100, RESEARCH/Physio Care Concept, , ch: SC, vol:0.5-10	yes	yes	yes	yes		
EPPENDORF@, 4710, DIGITAL, , ch: SC, vol:10-100					yes	
EPPENDORF@, 4710, DIGITAL, ULTRAMICRO, ch: SC, vol:0.5-10			yes	yes		
EPPENDORF@, 4810, ADJUSTABLE, , ch: SC, vol:10-100					yes	



Finnpipette Thermo Electron, , Digital, , ch: MC 12, vol:0.5-10	yes	yes				
Pipetter Description	5984D	5984C	5984F	5984H	5984L	5984N
Finnpipette Thermo Electron, , FP Focus Micro, , ch: SC, vol:3-30	yes	yes				
Finnpipette Thermo Electron, , FP Focus, , ch: MC 12, vol:1-10	yes	yes				
Finnpipette Thermo Electron, , FP Focus, , ch: SC, vol:0.3-3	yes	yes				
Finnpipette Thermo Electron, , Novus USA, , ch: MC 12 E, vol:1-10	yes	yes				
Finnpipette Thermo Electron, , Novus USA, , ch: MC 8 E, vol:1-10	yes	yes				
Finnpipette Thermo Scientific, , FP Digital, , ch: MC 8, vol:0.5-10	yes	yes				
Finnpipette Thermo Scientific, , FP Focus, , ch: MC 8, vol:1-10	yes	yes				
FINNPIPETTE®, , DIGITAL MCP'S, , ch: MC-8, vol:50-300					yes	
FINNPIPETTE®, , FOCUS THERMO LABSYSTEMS, , ch: MC-8, vol:1-10	yes	yes	yes	yes		
FINNPIPETTE®, , FOCUS THERMO LABSYSTEMS, , ch: MC-8, vol:5-50						
FINNPIPETTE®, , FOCUS THERMO LABSYSTEMS, , ch: SC, vol:100-1000						yes
FINNPIPETTE®, , FOCUS THERMO LABSYSTEMS, , ch: SC, vol:10-100					yes	
FINNPIPETTE®, , FOCUS THERMO LABSYSTEMS, , ch: SC, vol:1-10						
FINNPIPETTE®, , FOCUS THERMO LABSYSTEMS, , ch: SC, vol:30-300					yes	
FINNPIPETTE®, , FOCUS THERMO LABSYSTEMS, MICRO, ch: SC, vol:1-10	yes	yes				
FINNPIPETTE®, 4500, LABSYSTEMS, , ch: SC, vol:0.2-20						
FINNPIPETTE®, 4510, LABSYSTEMS, , ch: MC-12, vol:50-300						
FINNPIPETTE®, 4510, THERMO LABSYSTEMS, , ch: MC-12, vol:50-300						
FINNPIPETTE®, 4510, THERMO LABSYSTEMS, , ch: MC-8, vol:50-300						
FINNPIPETTE®, AS804, LABSYSTEMS DIGITAL MCP'S, , ch: MC-8, vol:5-50						
FISHERBRAND®, , FINNPIPETTE II, , ch: MC-8, vol:50-300						
FISHERBRAND®, , FINNPIPETTE II, , ch: SC, vol:0.2-2	yes	yes	yes	yes		
FISHERBRAND®, , FINNPIPETTE II, , ch: SC, vol:100-1000						yes
FISHERBRAND®, , FINNPIPETTE II, , ch: SC, vol:20-200					yes	
FISHERBRAND®, , FINNPIPETTE II, , ch: SC, vol:5-50					yes	



FISHERBRAND®, FINNPIPETTE II, MICRO, ch: SC, vol:0.5-10	yes	yes	yes	yes		
Pipetter Description	5984D	5984C	5984F	5984H	5984L	5984N
FISHERBRAND®, FINNPIPETTE, , ch: MC, vol:0.2-2	yes	yes	yes	yes		
FISHERBRAND®, FINNPIPETTE, , ch: MC-8, vol:0.5-10	yes	yes	yes	yes		
FISHERBRAND®, FINNPIPETTE, , ch: MC-8, vol:50-300					yes	
FISHERBRAND®, FINNPIPETTE, , ch: MC-8, vol:5-50						
FISHERBRAND®, FINNPIPETTE, , ch: SC, vol:0.2-2	yes	yes	yes	yes		
FISHERBRAND®, FINNPIPETTE, , ch: SC, vol:200-1000						yes
FISHERBRAND®, FINNPIPETTE, , ch: SC, vol:20-200					yes	
FISHERBRAND®, FINNPIPETTE, , ch: SC, vol:40-200					yes	
FISHERBRAND®, FINNPIPETTE, , ch: SC, vol:5-40					yes	
FISHERBRAND®, FINNPIPETTE, , ch: SC, vol:5-50					yes	
FISHERBRAND®, FINNPIPETTE, MICRO, ch: SC, vol:0.5-10	yes	yes	yes	yes		
GILSON, , Pipetman, M8, ch: MC8, vol:20-200						
GILSON, , Pipetman, P-10, ch: SC, vol:0.5-10	yes	yes				
GILSON, , Pipetman, P-100, ch: SC, vol:10-100						
GILSON, , Pipetman, P-1000, ch: SC, vol:100-1000						yes
GILSON, , Pipetman, P-2, ch: SC, vol:0.1-2	yes	yes				
GILSON, , Pipetman, P-20, ch: SC, vol:2-20			yes	yes		
GILSON, , Pipetman, P-200, ch: SC, vol:50-200					yes	
HAMILTON, , SOFTGRIP™, , ch: MC-8, vol:30-300					yes	
HAMILTON, , SOFTGRIP™, , ch: MC-8, vol:5-50					yes	
HAMILTON, , SOFTGRIP™, , ch: SC, vol:0.2-2	yes	yes	yes	yes		
HAMILTON, , SOFTGRIP™, , ch: SC, vol:100-1000						yes
HAMILTON, , SOFTGRIP™, , ch: SC, vol:10-100					yes	
HAMILTON, , SOFTGRIP™, , ch: SC, vol:1-10	yes	yes	yes	yes		
HAMILTON, , SOFTGRIP™, , ch: SC, vol:2.5-25					yes	



HAMILTON, , SOFTGRIP™, , ch: SC, vol:30-300						
Pipetter Description	5984D	5984C	5984F	5984H	5984L	5984N
ISC Bio Express, P-3947, GENEMATE®, , ch: SC, vol:0.1-2	yes	yes	yes	yes		
ISC Bio Express, P-3947, GENEMATE®, , ch: SC, vol:0.2-20						
ISC Bio Express, P-3947, GENEMATE®, , ch: SC, vol:0.5-10	yes	yes	yes	yes		
ISC Bio Express, P-3947, GENEMATE®, , ch: SC, vol:10-100					yes	
ISC Bio Express, P-3947, GENEMATE®, , ch: SC, vol:20-200					yes	
JENCONS, , , ACCURETTE, ch: SC, vol:1-10	yes	yes				
JENCONS, , , ACCURETTE, ch: SC, vol:20-100					yes	
JENCONS, , , ACCURETTE, ch: SC, vol:2-20						
JENCONS, , , ACCURETTE, ch: SC, vol:50-200					yes	
JENCONS, , , ACCURETTE, ch: SC, vol:50-200						
JENCONS, , SEALPETTE® ELECTRONIC, , ch: SC, E, vol:0.2-10	yes	yes				
JENCONS, , SEALPETTE® ELECTRONIC, , ch: SC, E, vol:5-100					yes	
JENCONS, , SEALPETTE®, , ch: MC-12, vol:50-250						
JENCONS, , SEALPETTE®, , ch: MC-8, vol:50-300						
JENCONS, , SEALPETTE®, , ch: SC, vol:50-200					yes	
JENCONS, , SEALPETTE®, , ch: SC, vol:5-50					yes	
LABNET, , LABPETTE DIGITAL, , ch: MC-8 & 12, vol:1-10		yes				
LABNET, , LABPETTE DIGITAL, , ch: MC-8 & 12, vol:20-200					yes	
LABNET, VE-10 LAB3, LABPETTE DIGITAL, VE-10, ch: SC, vol:0.5-10,	yes	yes				
LABNET, VE-100 LAB3, LABPETTE DIGITAL, VE-100, ch: SC, vol:10-100						
LABNET, VE-2 LAB3, LABPETTE DIGITAL, VE-2, ch: SC, vol:0.1-2	yes	yes				
LABNET, VE-20 LAB3, LABPETTE DIGITAL, VE-20, ch: SC, vol:0.2-20						
LABNET, VE-200 LAB3, LABPETTE DIGITAL, VE-200, ch: SC, vol:20-200					yes	
MATRIX, , MULTI-ELECTRAPETTE, , ch: MC, E, vol:0-25					yes	
METTLER-TOLEDO, , VOLUMATE, , ch: SC, vol:0.1-2.5	yes					



METTLER-TOLEDO, , VOLUMATE, , ch: SC, vol:0.5-10	yes	yes	yes	yes		
Pipetter Description	5984D	5984C	5984F	5984H	5984L	5984N
METTLER-TOLEDO, , VOLUMATE, , ch: SC, vol:10-100						
METTLER-TOLEDO, , VOLUMATE, , ch: SC, vol:20-200					yes	
METTLER-TOLEDO, , VOLUMATE, , ch: SC, vol:2-20						
NICHIRYO, , NICHIPET EX, , ch: SC, vol:10-100						
NICHIRYO, , NICHIPET EX, , ch: SC, vol:20-200						
NICHIRYO, , NICHIPET EX, , ch: SC, vol:2-20						
NICHIRYO, 5000, , , ch: SC, vol:10-50						
NICHIRYO, 5000, , , ch: SC, vol:200-1000						yes
NICHIRYO, 5000, , , ch: SC, vol:40-200					yes	
NICHIRYO, 5000F, FIXED VOLUME, , ch: SC, vol:10-10	yes	yes	yes	yes		
NICHIRYO, 5000F, FIXED VOLUME, , ch: SC, vol:200-200						
NICHIRYO, 5000V-L, , , ch: SC, vol:40-200						
NICHIRYO, 5000V-M, TRIPLE VOLUME, ch: SC, vol:200-1000200/500/1000µl						yes
NICHIRYO, 5000V-M, TRIPLE VOLUME, ch: SC, vol:50-20050/100/200µl					yes	
NICHIRYO, 7000, NICHIPET DIGITAL, , ch: MC-8, vol:40-200					yes	
NICHIRYO, 7000, NICHIPET DIGITAL, , ch: MC-8, vol:40-200						
OXFORD®, , BENCHMATE®, , ch: MC-12, vol:40-200					yes	
OXFORD®, , BENCHMATE®, , ch: SC, vol:0.5-10	yes	yes	yes	yes		
OXFORD®, , BENCHMATE®, , ch: SC, vol:10-50						
OXFORD®, , BENCHMATE®, , ch: SC, vol:200-1000						yes
OXFORD®, , BENCHMATE®, , ch: SC, vol:40-200					yes	
RAININ, , PIPET-PLUS, R-10, ch: SC, vol:0.5-10	yes	yes				
RAININ, , PIPET-PLUS, R-100, ch: SC, vol:10-100					yes	
RAININ, , PIPET-PLUS, R-1000, ch: SC, vol:100-1000						yes
RAININ, , PIPET-PLUS, R-2, ch: SC, vol:0.1-2	yes	yes				



RAININ, , PIPET-PLUS, R-20, ch: SC, vol:2-20			yes	yes		
Pipetter Description	5984D	5984C	5984F	5984H	5984L	5984N
RAININ, , PIPET-PLUS, R-200, ch: SC, vol:50-200					yes	
SOCOREX, 811, , , ch: SC, vol:100-100					yes	
SOCOREX, 811, , , ch: SC, vol:5-5			yes			
SOCOREX, 815, Acura Fixed Volume, , ch: SC, vol:200-200						
SOCOREX, 821, ch: SC, vol:200-1000						yes
SOCOREX, 821, , , ch: SC, vol:50-200					yes	
SOCOREX, 821, , , ch: SC, vol:5-50						
SOCOREX, 822, CALIBRA®, , ch: SC, vol:0.4-2			yes			
SOCOREX, 822, CALIBRA®, , ch: SC, vol:100-1000						yes
SOCOREX, 822, CALIBRA®, , ch: SC, vol:10-100						
SOCOREX, 822, CALIBRA®, , ch: SC, vol:1-10	yes	yes	yes	yes		
SOCOREX, 822, CALIBRA®, , ch: SC, vol:1-10			yes			
SOCOREX, 822, CALIBRA®, , ch: SC, vol:20-200					yes	
SOCOREX, 822, CALIBRA®, , ch: SC, vol:2-20						
SOCOREX, 825, Acura, , ch: SC, vol:0.1-2	yes	yes	yes	yes		
SOCOREX, 825, Acura, , ch: SC, vol:0.5-10	yes	yes	yes	yes		
SOCOREX, 825, Acura, , ch: SC, vol:100-1000						yes
SOCOREX, 825, Acura, , ch: SC, vol:10-100					yes	
SOCOREX, 825, Acura, , ch: SC, vol:20-200					yes	
SOCOREX, 825, Acura, , ch: SC, vol:2-20						
SOCOREX, 851, ACURA ANALOGIC, ch: MC-8, vol:50-200					yes	
SOCOREX, 851, ACURA ANALOGIC, , ch: MC-8,12, vol:5-50						
SOCOREX, 852, CALIBRA®, , ch: MC-8, vol:10-100						
SOCOREX, 852, CALIBRA®/WHEATON, , ch: MC-8, vol:1-10	yes	yes				
TITERTEK®/ PLUS ICN, , ch: MC-8, E, vol:5-100					yes	



TITERTEK® / PLUS, ch: MC-12, vol:5-50						
Pipetter Description	5984D	5984C	5984F	5984H	5984L	5984N
TITERTEK®, DIGITAL, ch: MC-12, vol:50-300					yes	
TITERTEK®, DIGITAL, ch: MC-12, vol:5-50						
TITERTEK®,DIGITAL, ch: SC, vol:200-1000						yes
TITERTEK®, DIGITAL, ch: SC, vol:40-200					yes	
TITERTEK®,DIGITAL, ch: SC, vol:5-40					yes	
VISTALAB,MLA®, ch: SC, vol:100-100, 1025					yes	
VISTALAB, , MLA®, , ch: SC, vol:20-20					yes	
VISTALAB, , MLA®, , ch: SC, vol:20-20						
VISTALAB, , MLA®, , ch: SC, vol:50-200						
VISTALAB, , MLA®, , ch: SC, vol:50-50					yes	
VISTALAB,MLA®, Precision D Calibratable, ch: SC, vol:10-10					yes	
VISTALAB,MLA®, Precision D Calibratable, ch:SC, vol:200-200						
VISTALAB,MLA®, Precision D Calibratable, ch:SC, vol:25-25					yes	
VISTALAB,MLA®, Precision D Calibratable, ch: SC, vol:5-5					yes	
VISTALAB,OVATION, ch: MC-12, vol:5-250						
VISTALAB, OVATION, ch: SC, vol:0.2-2	yes					
VISTALAB,OVATION,ch: SC, vol:10-100					yes	
VISTALAB,OVATION,ch: SC, vol:1-10		yes	yes	yes		
VISTALAB,OVATION,ch: SC, vol:2-20		yes	yes	yes		
VISTALAB, , OVATION, gray cone installed, ch: SC, vol:20-200						
VISTALAB, , OVATION, OVATION, ch: MC-8, vol:5-250						
VISTALAB, , OVATION, white cone installed, ch: SC, vol:100-1000						yes
VISTALAB, , OVATION-white cone, old version, white cone, ch: SC, vol:20-200						
VWR Signature, , International, , ch: MC 12 E, vol:10-200						
VWR Signature, , International, , ch: MC 12, vol:50-300						



VWR Signature, , International, , ch: MC 12, vol:5-50						
Pipetter Description	5984D	5984C	5984F	5984H	5984L	5984N
VWR Signature, , International, , ch: MC 8, vol:1-10	yes	yes	yes	yes		
VWR Signature, , International, , ch: SC, vol:0.1-2	yes	yes	yes	yes		
VWR Signature, , International, , ch: SC, vol:0.5-10	yes	yes				
VWR Signature, , International, , ch: SC, vol:10-100					yes	
VWR Signature, , International, , ch: SC, vol:20-200					yes	
VWR Signature, , International, , ch: SC, vol:2-20						