

Kompatibilitätsübersicht der Fisherbrand SureOne™ Spitzen

Diese Liste erhebt keinen Anspruch auf Vollständigkeit. Sollten Sie Ihre Pipette in dieser Liste nicht finden, kontaktieren Sie uns bitte unter der Tel. 07851-99457-0 oder per e-mail : sales@novodirect.de

Pipetter Description	5982C 5982D 5982E	59824 59826 59829 5982B	5982F 5982H 5982K	59835 59837 59838 59839 5983A 5983C 5983D	5982L 5982N 5982T 5982Y 59831 59833 59834	5983E 5983F 5983H 5983K 5983L 5983M 5983N	5983T 5983W 5983X 59840	59841 59842 59843 59844 59845 59846 59847 59848
BIOHIT, , eLine, e10, ch: MC8, vol:0.2-10	yes	yes						
BIOHIT, , eLine, e10, ch: SC, E, vol:0.2-10	yes	yes						
BIOHIT, , eLine, e120, ch: MC8, vol:5-120					yes		yes	
BIOHIT, , eLine, e120, ch: SC, E, vol:5-120					yes			
BIOHIT, , mLINE®, m10, ch: MC12, vol:0.5-10	yes	yes						
BIOHIT, , mLINE®, m10, ch: MC8, vol:0.5-10	yes	yes						
BIOHIT, , mLINE®, m10, ch: SC, vol:0.5-10	yes	yes						
BIOHIT, , mLINE®, m100, ch: MC12, vol:5-100					yes		yes	
BIOHIT, , mLINE®, m100, ch: MC8, vol:5-100					yes		yes	
BIOHIT, , mLINE®, m100, ch: SC, vol:10-100					yes		yes	
BIOHIT, , mLINE®, m1000, ch: SC, vol:100-1000								yes
BIOHIT, , mLINE®, m20, ch: SC, vol:2-20					yes			



BIOHIT, , mLINE®, m200, ch: SC, vol:20-200					yes		yes	
				59835	5982L	5983E		59841
				59837	5982N	5983F		59842
				59838	5982T	5983H		59843
		59824		59839	5982Y	5983K	5983T	59844
	5982C	59826	5982F	5983A	59831	5983L	5983W	59845
	5982D	59829	5982H	5983C	59833	5983M	5983X	59846
	5982E	5982B	5982K	5983D	59834	5983N	59840	59847
								59848
BIOHIT, , mLINE®, m300, ch: MC8, vol:30-300					yes			
BIOHIT, , PROLINE® ELECTRONIC, , ch: MC-8, E, vol:25-250				yes	yes	yes		
BIOHIT, , PROLINE® ELECTRONIC, , ch: MC-8, E, vol:5-100				yes	yes	yes		
BIOHIT, , PROLINE® ELECTRONIC, , ch: SC, E, vol:10-250							yes	
BIOHIT, , PROLINE® ELECTRONIC, , ch: SC, E, vol:5-100				yes	yes	yes		
BIOHIT, , PROLINE®, , ch: MC, vol:5-50					yes			
BIOHIT, , PROLINE®, , ch: MC-12, vol:50-250							yes	
BIOHIT, , PROLINE®, , ch: MC-8, vol:50-300							yes	
BIOHIT, , PROLINE®, , ch: SC, vol:0.5-10				yes	yes	yes	yes	
BIOHIT, , PROLINE®, , ch: SC, vol:5-50				yes	yes	yes	yes	
BIOHIT, , PROLINE®, PROLINE®, ch: SC, vol:100-1000								yes
BIOHIT, , PROLINE®, PROLINE®, ch: SC, vol:20-200				yes	yes	yes	yes	
BIOZYM, , PRECISION, , ch: SC, E, vol:0.2-10	yes	yes						
BRANDTECH, , TRANSFERPETTE®, , ch: MC-12, vol:10-100				yes	yes	yes		
BRANDTECH, , TRANSFERPETTE®, , ch: MC-12, vol:2.5-25				yes		yes		
BRANDTECH, , TRANSFERPETTE®, , ch: MC-12, vol:5-50				yes	yes	yes		
BRANDTECH, , TRANSFERPETTE®, , ch: SC, vol:0.5-10	yes	yes	yes					
BRANDTECH, , TRANSFERPETTE®, , ch: SC, vol:10-100				yes	yes	yes		
BRANDTECH, , TRANSFERPETTE®, , ch: SC, vol:20-200				yes	yes	yes	yes	
BRANDTECH, , TRANSFERPETTE®, , ch: SC, vol:5-50				yes	yes	yes		



BRINKMANN, , TRANSFERPETTE®, , ch: MC-12, vol:5-25				yes		yes		
				59835	5982L	5983E		59841
				59837	5982N	5983F		59842
				59838	5982T	5983H		59843
		59824		59839	5982Y	5983K	5983T	59844
	5982C	59826	5982F	5983A	59831	5983L	5983W	59845
	5982D	59829	5982H	5983C	59833	5983M	5983X	59846
	5982E	5982B	5982K	5983D	59834	5983N	59840	59847
								59848
BRINKMANN, , TRANSFERPETTE®, , ch: MC-8, vol:50-200					yes			
CappAero, , , CappAero, ch: MC-8, vol:25-200					yes			
CappAero, , , CappAero, ch: SC, vol:0.5-10	yes	yes						
CappAero, , , CappAero, ch: SC, vol:25-200					yes			
CappAero, , , Softline, ch: MC8, vol:5-50					yes			
CappAero, , , Softline, ch: SC, vol:0.5-10	yes	yes						
CappAero, , , Softline, ch: SC, vol:10-100				yes	yes	yes		
CappAero, , , Softline, ch: SC, vol:25-200				yes	yes	yes		
CappAero, , , Softline, ch: SC, vol:5-50				yes	yes	yes		
COSTAR®, , 8 PETTE, , ch: MC8, vol:20-200				yes	yes	yes	yes	
COSTAR®, , ONEPETTE™, , ch: SC, vol:10-100				yes	yes	yes	yes	
COSTAR®, , ONEPETTE™, , ch: SC, vol:1-20				yes	yes	yes	yes	
COSTAR®, , ONEPETTE™, , ch: SC, vol:20-200				yes	yes	yes	yes	
EPPENDORF®, , REFERENCE, Physio Care, ch: SC, vol:100-1000								yes
EPPENDORF®, , REFERENCE, Physio Care, ch: SC, vol:50-200				yes	yes	yes	yes	
EPPENDORF®, , REFERENCE, Physio Care/Metal Cone, ch: SC, vol:2-20	yes	yes	yes					
EPPENDORF®, , REFERENCE, Physio Care/Metal Cone, ch: SC, vol:5-10	yes	yes	yes					
EPPENDORF®, , Research pro, Physio Care, ch: MC-8, vol:50-1200								yes
EPPENDORF®, , RESEARCH, Physio Care, ch: MC 12, vol:10-100				yes		yes		
EPPENDORF®, , RESEARCH, Physio Care, ch: SC, vol:0.1-2.5	yes	yes						



GILSON, , Pipetman, P-100, ch: SC, vol:10-100				yes	yes	yes	yes	
				59835	5982L	5983E		59841
				59837	5982N	5983F		59842
				59838	5982T	5983H		59843
				59839	5982Y	5983K		59844
				5983A	59831	5983L	5983T	59845
				5983C	59833	5983M	5983W	59846
				5983D	59834	5983N	5983X	59847
							59840	59848
Pipetter Description	5982C	59824	5982F					
	5982D	59826	5982H					
	5982E	59829	5982K					
		5982B						
GILSON, , Pipetman, P-1000, ch: SC, vol:100-1000								yes
GILSON, , Pipetman, P-2, ch: SC, vol:0.1-2	yes	yes						
GILSON, , Pipetman, P-20, ch: SC, vol:2-20			yes	yes	yes	yes		
GILSON, , Pipetman, P-200, ch: SC, vol:50-200				yes	yes	yes	yes	
HAMILTON, , SOFTGRIP™, , ch: MC-8, vol:30-300				yes		yes	yes	
HAMILTON, , SOFTGRIP™, , ch: MC-8, vol:5-50				yes	yes	yes	yes	
HAMILTON, , SOFTGRIP™, , ch: SC, vol:0.2-2	yes	yes	yes					
HAMILTON, , SOFTGRIP™, , ch: SC, vol:100-1000								yes
HAMILTON, , SOFTGRIP™, , ch: SC, vol:10-100				yes	yes	yes	yes	
HAMILTON, , SOFTGRIP™, , ch: SC, vol:1-10	yes	yes	yes					
HAMILTON, , SOFTGRIP™, , ch: SC, vol:2.5-25				yes	yes	yes	yes	
HAMILTON, , SOFTGRIP™, , ch: SC, vol:30-300					yes		yes	
ISC Bio Express, P-3947, GENEMATE®, , ch: SC, vol:0.1-2	yes	yes	yes					
ISC Bio Express, P-3947, GENEMATE®, , ch: SC, vol:0.2-20					yes	yes		
ISC Bio Express, P-3947, GENEMATE®, , ch: SC, vol:0.5-10	yes	yes	yes					
ISC Bio Express, P-3947, GENEMATE®, , ch: SC, vol:10-100				yes	yes	yes	yes	
ISC Bio Express, P-3947, GENEMATE®, , ch: SC, vol:20-200				yes	yes	yes	yes	
JENCONS, , , ACCURETTE, ch: SC, vol:100-1000								yes
JENCONS, , , ACCURETTE, ch: SC, vol:1-10	yes	yes						
JENCONS, , , ACCURETTE, ch: SC, vol:20-100				yes	yes	yes	yes	



JENCONS, , , ACCURETTE, ch: SC, vol:2-20				yes		yes		
				59835	5982L	5983E		59841
				59837	5982N	5983F		59842
				59838	5982T	5983H		59844
				59839	5982Y	5983K	5983T	59845
	5982C	59824	5982F	5983A	59831	5983L	5983W	59846
	5982D	59826	5982H	5983C	59833	5983M	5983X	59847
	5982E	5982B	5982K	5983D	59834	5983N	59840	59848
JENCONS, , , ACCURETTE, ch: SC, vol:50-200				yes	yes	yes	yes	
JENCONS, , SEALPETTE® ELECTRONIC, , ch: SC, E, vol:0.2-10	yes	yes						
JENCONS, , SEALPETTE® ELECTRONIC, , ch: SC, E, vol:10-250					yes		yes	
JENCONS, , SEALPETTE® ELECTRONIC, , ch: SC, E, vol:5-100				yes	yes	yes	yes	
JENCONS, , SEALPETTE®, , ch: MC-12, vol:50-250					yes		yes	
JENCONS, , SEALPETTE®, , ch: MC-8, vol:50-300					yes		yes	
JENCONS, , SEALPETTE®, , ch: SC, vol:0.5-10				yes	yes	yes	yes	
JENCONS, , SEALPETTE®, , ch: SC, vol:50-200				yes	yes	yes	yes	
JENCONS, , SEALPETTE®, , ch: SC, vol:5-50				yes	yes	yes	yes	
LABNET, , LABPETTE DIGITAL, , ch: MC-8 & 12, vol:1-10		yes						
LABNET, , LABPETTE DIGITAL, , ch: MC-8 & 12, vol:20-200				yes	yes	yes		
LABNET, VE-10 LAB3, LABPETTE DIGITAL, VE-10, ch: SC, vol:0.5-10	yes	yes						
LABNET, VE-100 LAB3, LABPETTE DIGITAL, VE-100, ch: SC, vol:10-100				yes	yes	yes		
LABNET, VE-2 LAB3, LABPETTE DIGITAL, VE-2, ch: SC, vol:0.1-2	yes	yes						
LABNET, VE-20 LAB3, LABPETTE DIGITAL, VE-20, ch: SC, vol:0.2-20				yes	yes	yes		
LABNET, VE-200 LAB3, LABPETTE DIGITAL, VE-200, ch: SC, vol:20-200				yes		yes		
MATRIX, , MULTI-ELECTRAPETTE, , ch: MC, E, vol:0-25					yes		yes	
MATRIX, , MULTI-ELECTRAPETTE, , ch: MC-8, E, vol:5-250					yes		yes	
METTLER-TOLEDO, , LIQUIDSYSTEMS, , ch: MC-8, vol:5-50					yes			
METTLER-TOLEDO, , LIQUIDSYSTEMS, , ch: SC, vol:10-100					yes			



SOCOREX, 821, , , ch: SC, vol:50-200				yes		yes		
				59835	5982L	5983E		59841
				59837	5982N	5983F		59842
				59838	5982T	5983H		59844
				59839	5982Y	5983K	5983T	59845
	5982C	59824	5982F	5983A	59831	5983L	5983W	59846
	5982D	59826	5982H	5983C	59833	5983M	5983X	59847
	5982E	59829	5982K	5983D	59834	5983N	59840	59848
		5982B						
SOCOREX, 821, , , ch: SC, vol:5-50				yes		yes		
SOCOREX, 822, CALIBRA®, , ch: SC, vol:0.4-2			yes					
SOCOREX, 822, CALIBRA®, , ch: SC, vol:100-1000								yes
SOCOREX, 822, CALIBRA®, , ch: SC, vol:10-100				yes	yes	yes		
SOCOREX, 822, CALIBRA®, , ch: SC, vol:1-10	yes	yes	yes					
SOCOREX, 822, CALIBRA®, , ch: SC, vol:1-10			yes					
SOCOREX, 822, CALIBRA®, , ch: SC, vol:20-200				yes		yes		
SOCOREX, 822, CALIBRA®, , ch: SC, vol:2-20				yes		yes		
SOCOREX, 825, Acura, , ch: SC, vol:0.1-2	yes	yes	yes					
SOCOREX, 825, Acura, , ch: SC, vol:0.5-10	yes	yes	yes					
SOCOREX, 825, Acura, , ch: SC, vol:100-1000								yes
SOCOREX, 825, Acura, , ch: SC, vol:10-100				yes	yes	yes	yes	
SOCOREX, 825, Acura, , ch: SC, vol:20-200				yes	yes	yes	yes	
SOCOREX, 825, Acura, , ch: SC, vol:2-20				yes	yes	yes	yes	
SOCOREX, 851, ACURA ANALOGIC, , ch: MC-8, vol:50-200				yes	yes	yes		
SOCOREX, 851, ACURA ANALOGIC, , ch: MC-8,12, vol:5-50				yes		yes		
SOCOREX, 852, CALIBRA®, , ch: MC-8, vol:10-100,				yes	yes	yes		
SOCOREX, 852, CALIBRA®/WHEATON, , ch: MC-8, vol:1-10	yes	yes						
SOCOREX, 855, Acura, , ch: MC12, vol:20-200					yes			
SOCOREX, 855, Acura, , ch: MC12, vol:40-350					yes			



TITERTEK®, , , ch: MC-12, vol:50-200				yes	yes	yes	yes	
				59835	5982L	5983E		59841
				59837	5982N	5983F		59842
				59838	5982T	5983H		59844
				59839	5982Y	5983K	5983T	59845
	5982C	59824	5982F	5983A	59831	5983L	5983W	59846
	5982D	59826	5982H	5983C	59833	5983M	5983X	59847
	5982E	59829	5982K	5983D	59834	5983N	59840	59848
TITERTEK®, , , ch: MC-8, vol:50-300				yes	yes	yes	yes	
TITERTEK®, , , ch: MC-8, vol:5-50							yes	
TITERTEK®, , , ch: SC, vol:0.5-10				yes	yes			
TITERTEK®, / PLUS ICN, , ch: MC-8, E, vol:5-100				yes	yes	yes		
TITERTEK®, / PLUS, , ch: MC-12, vol:5-50				yes		yes		
TITERTEK®, / PLUS, , ch: SC, vol:100-100				yes	yes	yes		
TITERTEK®, / PLUS, , ch: SC, vol:250-250								yes
TITERTEK®, / PLUS, , ch: SC, vol:5-5			yes					
TITERTEK®, DIGITAL, , ch: MC-12, vol:50-300				yes	yes	yes	yes	
TITERTEK®, DIGITAL, , ch: MC-12, vol:5-50				yes		yes	yes	
TITERTEK®, DIGITAL, , ch: SC, vol:0.5-10				yes	yes	yes	yes	
TITERTEK®, DIGITAL, , ch: SC, vol:200-1000								yes
TITERTEK®, DIGITAL, , ch: SC, vol:40-200				yes	yes	yes	yes	
TITERTEK®, DIGITAL, , ch: SC, vol:5-40				yes	yes	yes	yes	
VISTALAB, , MLA®, , ch: SC, vol:100-100, 1025				yes	yes	yes	yes	
VISTALAB, , MLA®, , ch: SC, vol:20-20				yes	yes	yes	yes	
VISTALAB, , MLA®, , ch: SC, vol:50-200							yes	
VISTALAB, , MLA®, , ch: SC, vol:50-50				yes	yes	yes	yes	
VISTALAB, , MLA®, Precision D Calibratable, ch: SC, vol:10-10				yes	yes	yes	yes	
VISTALAB, , MLA®, Precision D Calibratable, ch: SC, vol:200-200				yes	yes	yes	yes	



VISTALAB, , MLA®, Precision D Calibratable, ch: SC, vol:25-25				yes	yes	yes	yes	
				59835	5982L	5983E		59841
				59837	5982N	5983F		59842
				59838	5982T	5983H		59844
				59839	5982Y	5983K	5983T	59845
	5982C	59824	5982F	5983A	59831	5983L	5983W	59846
	5982D	59826	5982H	5983C	59833	5983M	5983X	59847
	5982E	5982B	5982K	5983D	59834	5983N	59840	59848
VISTALAB, , MLA®, Precision D Calibratable, ch: SC, vol:5-5				yes	yes	yes	yes	
VISTALAB, , OVATION, , ch: MC-12, vol:5-250					yes			
VISTALAB, , OVATION, , ch: SC, vol:0.2-2,	yes	yes						
VISTALAB, , OVATION, , ch: SC, vol:10-100				yes	yes	yes	yes	
VISTALAB, , OVATION, , ch: SC, vol:1-10		yes	yes					
VISTALAB, , OVATION, , ch: SC, vol:2-20		yes	yes					
VISTALAB, , OVATION, gray cone installed, ch: SC, vol:20-200					yes			
VISTALAB, , OVATION, OVATION, ch: MC-8, vol:5-250					yes			
VISTALAB, , OVATION, white cone installed, ch: SC, vol:100-1000								yes
VISTALAB,OVATION-white cone, old version ch:SC, vol:20-200					yes		yes	
VWR Signature, , International, , ch: MC 12 E, vol:10-200					yes			
VWR Signature, , International, , ch: MC 12, vol:50-300					yes			
VWR Signature, , International, , ch: MC 12, vol:5-50				yes	yes	yes		
VWR Signature, , International, , ch: MC 8, vol:1-10	yes	yes	yes					
VWR Signature, , International, , ch: SC, vol:0.1-2	yes	yes	yes					
VWR Signature, , International, , ch: SC, vol:0.5-10	yes	yes						
VWR Signature, , International, , ch: SC, vol:100-1000								yes
VWR Signature, , International, , ch: SC, vol:10-100				yes	yes	yes	yes	
VWR Signature, , International, , ch: SC, vol:20-200				yes	yes	yes	yes	
VWR Signature, , International, , ch: SC, vol:2-20				yes	yes	yes	yes	



Fisherbrand[®]

Tableau de compatibilité pointes Fisherbrand SureOne Avec Filtre

Ci joint un tableau comparatif qui a été testé avec les pointes fisherbrand SureOne. Dans le cas où la micropipette du client n'a pas été testée avec les pointes SureOne, n'hésitez pas à faire une demande d'échantillon de pointes auprès de MM ou FC.

Pipetter Description	5984D	5984C	5984F	5984H	5984L	5984N
BIOHIT, , eLine, e10, ch: MC8, vol:0.2-10	yes	yes				
BIOHIT, , eLine, e10, ch: SC, E, vol:0.2-10	yes	yes				
BIOHIT, , eLine, e300, ch: MC8 E, vol:10-300						
BIOHIT, , mLINE®, m10, ch: MC12, vol:0.5-10	yes	yes				
BIOHIT, , mLINE®, m10, ch: MC8, vol:0.5-10	yes	yes				
BIOHIT, , mLINE®, m10, ch: SC, vol:0.5-10	yes	yes				
BIOHIT, , mLINE®, m100, ch: SC, vol:10-100					yes	
BIOHIT, , mLINE®, m1000, ch: SC, vol:100-1000						yes
BIOHIT, , mLINE®, m200, ch: SC, vol:20-200						
BIOHIT, , mLINE®, m300, ch: MC12, vol:30-300						
BIOHIT, , mLINE®, m300, ch: MC8, vol:30-300						
BIOHIT, , PROLINE® ELECTRONIC, , ch: MC-8, E, vol:25-250					yes	
BIOHIT, , PROLINE® ELECTRONIC, , ch: MC-8, E, vol:5-100					yes	
BIOHIT, , PROLINE® ELECTRONIC, , ch: SC, E, vol:5-100					yes	
BIOHIT, , PROLINE®, , ch: SC, vol:5-50					yes	
BIOHIT, , PROLINE®, PROLINE®, ch: SC, vol:100-1000						yes



BIOHIT, , PROLINE®, PROLINE®, ch: SC, vol:20-200					yes	
Pipetter Description	5984D	5984C	5984F	5984H	5984L	5984N
BIOZYM, , PRECISION, , ch: SC, E, vol:0.2-10	yes	yes				
BRANDTECH, , TRANSFERPETTE®, , ch: MC-12, vol:10-100						
BRANDTECH, , TRANSFERPETTE®, , ch: MC-12, vol:5-50						
BRANDTECH, , TRANSFERPETTE®, , ch: SC, vol:0.5-10	yes	yes	yes	yes		
BRANDTECH, , TRANSFERPETTE®, , ch: SC, vol:10-100						
BRANDTECH, , TRANSFERPETTE®, , ch: SC, vol:20-200					yes	
BRANDTECH, , TRANSFERPETTE®, , ch: SC, vol:5-50						
BRINKMANN, , TRANSFERPETTE®, , ch: MC-12, vol:5-25						
BRINKMANN, , TRANSFERPETTE®, , ch: MC-8, vol:50-200					yes	
CappAero, , , CappAero, ch: MC-8, vol:25-200						
CappAero, , , CappAero, ch: SC, vol:0.5-10	yes	yes				
CappAero, , , CappAero, ch: SC, vol:25-200					yes	
CappAero, , , Softline, ch: MC8, vol:50-300						
CappAero, , , Softline, ch: MC8, vol:5-50						
CappAero, , , Softline, ch: SC, vol:0.5-10	yes	yes				
CappAero, , , Softline, ch: SC, vol:10-100					yes	
CappAero, , , Softline, ch: SC, vol:25-200					yes	
CappAero, , , Softline, ch: SC, vol:5-50					yes	
CappAero, , , Softline, ch: SC, vol:5-50						
COSTAR®, , 8 PETTE, , ch: MC8, vol:20-200					yes	
COSTAR®, , ONEPETTE™, , ch: SC, vol:10-100						
COSTAR®, , ONEPETTE™, , ch: SC, vol:1-20						
COSTAR®, , ONEPETTE™, , ch: SC, vol:20-200					yes	
EPENDORF®,REFERENCE, Physio Care, ch:SC,vol:100-1000						yes



EPPENDORF®,REFERENCE, Physio Care, ch: SC,vol:50-200					yes	
EPPENDORF®,REFERENCE, Physio Care/Metal Cone,ch:SC,vol:2-20	yes	yes	yes	yes		
Pipetter Description	5984D	5984C	5984F	5984H	5984L	5984N
EPPENDORF®,REFERENCE,Physio Care/Metal Cone, ch:SC, vol:5-10	yes	yes	yes	yes		
EPPENDORF®,RESEARCH, Physio Care, ch: MC 12, vol:10-100					yes	
EPPENDORF®,RESEARCH, Physio Care, ch: SC, vol:0.1-2.5	yes					
EPPENDORF®,RESEARCH, Physio Care, ch:SC,vol:100-1000						yes
EPPENDORF®,RESEARCH, Physio Care, ch:SC,vol:10-100					yes	
EPPENDORF®,RESEARCH,Physio Care, ch:SC, vol:20-200					yes	
EPPENDORF®, , RESEARCH, Physio Care, ch: SC, vol:2-20					yes	
EPPENDORF®, 2000, REFERENCE, , ch: SC, vol:0.1-2.5	yes					
EPPENDORF®, 2000, REFERENCE, , ch: SC, vol:0.5-10	yes	yes	yes	yes		
EPPENDORF®, 2000, REFERENCE, , ch: SC, vol:10-100					yes	
EPPENDORF®, 2000, REFERENCE, ch: SC, vol:200-1000						yes
EPPENDORF®, 2000, REFERENCE, , ch: SC, vol:50-200					yes	
EPPENDORF®, 2000, REFERENCE, GRAY, ch: SC, vol:2-20	yes	yes	yes	yes		
EPPENDORF®, 2000, REFERENCE, YELLOW, ch: SC, vol:2-20					yes	
EPPENDORF®, 2100, RESEARCH, , ch: MC-8, vol:0.5-10			yes			
EPPENDORF®, 2100, RESEARCH, , ch: MC-8, vol:10-100					yes	
EPPENDORF®, 2100, RESEARCH, , ch: SC, vol:0.1-2.5	yes					
EPPENDORF®, 2100, RESEARCH, , ch: SC, vol:100-1000						yes
EPPENDORF®, 2100, RESEARCH, , ch: SC, vol:10-100					yes	
EPPENDORF®, 2100, RESEARCH, , ch: SC, vol:20-200					yes	
EPPENDORF®, 2100, RESEARCH, , ch: SC, vol:2-20						
EPPENDORF®, 2100, RESEARCH/Physio Care Concept, , ch: SC, vol:0.5-10	yes	yes	yes	yes		
EPPENDORF®, 4710, DIGITAL, , ch: SC, vol:10-100					yes	
EPPENDORF®, 4710, DIGITAL, ULTRAMICRO, ch: SC, vol:0.5-10			yes	yes		



EPPENDORF®, 4810, ADJUSTABLE, , ch: SC, vol:10-100					yes	
Finnpipette Thermo Electron, , Digital, , ch: MC 12, vol:0.5-10	yes	yes				
Pipetter Description	5984D	5984C	5984F	5984H	5984L	5984N
Finnpipette Thermo Electron, , FP Focus Micro, , ch: SC, vol:3-30	yes	yes				
Finnpipette Thermo Electron, , FP Focus, , ch: MC 12, vol:1-10	yes	yes				
Finnpipette Thermo Electron, , FP Focus, , ch: SC, vol:0.3-3	yes	yes				
Finnpipette Thermo Electron, , Novus USA, , ch: MC 12 E, vol:1-10	yes	yes				
Finnpipette Thermo Electron, , Novus USA, , ch: MC 8 E, vol:1-10	yes	yes				
Finnpipette Thermo Scientific, , FP Digital, , ch: MC 8, vol:0.5-10	yes	yes				
Finnpipette Thermo Scientific, , FP Focus, , ch: MC 8, vol:1-10	yes	yes				
FINNPIPETTE®, , DGITAL MCP'S, , ch: MC-8, vol:50-300					yes	
FINNPIPETTE®, , FOCUS THERMO LABSYSTEMS, , ch: MC-8, vol:1-10	yes	yes	yes	yes		
FINNPIPETTE®, , FOCUS THERMO LABSYSTEMS, , ch: MC-8, vol:5-50						
FINNPIPETTE®, , FOCUS THERMO LABSYSTEMS, , ch: SC, vol:100-1000						yes
FINNPIPETTE®, , FOCUS THERMO LABSYSTEMS, , ch: SC, vol:10-100					yes	
FINNPIPETTE®, , FOCUS THERMO LABSYSTEMS, , ch: SC, vol:1-10						
FINNPIPETTE®, , FOCUS THERMO LABSYSTEMS, , ch: SC, vol:30-300					yes	
FINNPIPETTE®, , FOCUS THERMO LABSYSTEMS, MICRO, ch: SC, vol:1-10	yes	yes				
FINNPIPETTE®, 4500, LABSYSTEMS, , ch: SC, vol:0.2-20						
FINNPIPETTE®, 4510, LABSYSTEMS, , ch: MC-12, vol:50-300						
FINNPIPETTE®, 4510, THERMO LABSYSTEMS, , ch: MC-12, vol:50-300						
FINNPIPETTE®, 4510, THERMO LABSYSTEMS, , ch: MC-8, vol:50-300						
FINNPIPETTE®, AS804, LABSYSTEMS DGITAL MCP'S, , ch: MC-8, vol:5-50						
FISHERBRAND®, , FINNPIPETTE II, , ch: MC-8, vol:50-300						
FISHERBRAND®, , FINNPIPETTE II, , ch: SC, vol:0.2-2	yes	yes	yes	yes		
FISHERBRAND®, , FINNPIPETTE II, , ch: SC, vol:100-1000						yes
FISHERBRAND®, , FINNPIPETTE II, , ch: SC, vol:20-200					yes	



FISHERBRAND®, FINNPIPETTE II, , ch: SC, vol:5-50					yes	
FISHERBRAND®, FINNPIPETTE II, MICRO, ch: SC, vol:0.5-10	yes	yes	yes	yes		
Pipetter Description	5984D	5984C	5984F	5984H	5984L	5984N
FISHERBRAND®, FINNPIPETTE, , ch: MC, vol:0.2-2	yes	yes	yes	yes		
FISHERBRAND®, FINNPIPETTE, , ch: MC-8, vol:0.5-10	yes	yes	yes	yes		
FISHERBRAND®, FINNPIPETTE, , ch: MC-8, vol:50-300					yes	
FISHERBRAND®, FINNPIPETTE, , ch: MC-8, vol:5-50						
FISHERBRAND®, FINNPIPETTE, , ch: SC, vol:0.2-2	yes	yes	yes	yes		
FISHERBRAND®, FINNPIPETTE, , ch: SC, vol:200-1000						yes
FISHERBRAND®, FINNPIPETTE, , ch: SC, vol:20-200					yes	
FISHERBRAND®, FINNPIPETTE, , ch: SC, vol:40-200					yes	
FISHERBRAND®, FINNPIPETTE, , ch: SC, vol:5-40					yes	
FISHERBRAND®, FINNPIPETTE, , ch: SC, vol:5-50					yes	
FISHERBRAND®, FINNPIPETTE, MICRO, ch: SC, vol:0.5-10	yes	yes	yes	yes		
GILSON, , Pipetman, M8, ch: MC8, vol:20-200						
GILSON, , Pipetman, P-10, ch: SC, vol:0.5-10	yes	yes				
GILSON, , Pipetman, P-100, ch: SC, vol:10-100						
GILSON, , Pipetman, P-1000, ch: SC, vol:100-1000						yes
GILSON, , Pipetman, P-2, ch: SC, vol:0.1-2	yes	yes				
GILSON, , Pipetman, P-20, ch: SC, vol:2-20			yes	yes		
GILSON, , Pipetman, P-200, ch: SC, vol:50-200					yes	
HAMILTON, , SOFTGRIP™, , ch: MC-8, vol:30-300					yes	
HAMILTON, , SOFTGRIP™, , ch: MC-8, vol:5-50					yes	
HAMILTON, , SOFTGRIP™, , ch: SC, vol:0.2-2	yes	yes	yes	yes		
HAMILTON, , SOFTGRIP™, , ch: SC, vol:100-1000						yes
HAMILTON, , SOFTGRIP™, , ch: SC, vol:10-100					yes	
HAMILTON, , SOFTGRIP™, , ch: SC, vol:1-10	yes	yes	yes	yes		



HAMILTON, , SOFTGRIP™, , ch: SC, vol:2.5-25					yes	
HAMILTON, , SOFTGRIP™, , ch: SC, vol:30-300						
Pipetter Description	5984D	5984C	5984F	5984H	5984L	5984N
ISC Bio Express, P-3947, GENEMATE®, , ch: SC, vol:0.1-2	yes	yes	yes	yes		
ISC Bio Express, P-3947, GENEMATE®, , ch: SC, vol:0.2-20						
ISC Bio Express, P-3947, GENEMATE®, , ch: SC, vol:0.5-10	yes	yes	yes	yes		
ISC Bio Express, P-3947, GENEMATE®, , ch: SC, vol:10-100					yes	
ISC Bio Express, P-3947, GENEMATE®, , ch: SC, vol:20-200					yes	
JENCONS, , , ACCURETTE, ch: SC, vol:1-10	yes	yes				
JENCONS, , , ACCURETTE, ch: SC, vol:20-100					yes	
JENCONS, , , ACCURETTE, ch: SC, vol:2-20						
JENCONS, , , ACCURETTE, ch: SC, vol:50-200					yes	
JENCONS, , , ACCURETTE, ch: SC, vol:50-200						
JENCONS, , SEALPETTE® ELECTRONIC, , ch: SC, E, vol:0.2-10	yes	yes				
JENCONS, , SEALPETTE® ELECTRONIC, , ch: SC, E, vol:5-100					yes	
JENCONS, , SEALPETTE®, , ch: MC-12, vol:50-250						
JENCONS, , SEALPETTE®, , ch: MC-8, vol:50-300						
JENCONS, , SEALPETTE®, , ch: SC, vol:50-200					yes	
JENCONS, , SEALPETTE®, , ch: SC, vol:5-50					yes	
LABNET, , LABPETTE DIGITAL, , ch: MC-8 & 12, vol:1-10		yes				
LABNET, , LABPETTE DIGITAL, , ch: MC-8 & 12, vol:20-200					yes	
LABNET, VE-10 LAB3, LABPETTE DIGITAL, VE-10, ch: SC, vol:0.5-10,	yes	yes				
LABNET, VE-100 LAB3, LABPETTE DIGITAL, VE-100, ch: SC, vol:10-100						
LABNET, VE-2 LAB3, LABPETTE DIGITAL, VE-2, ch: SC, vol:0.1-2	yes	yes				
LABNET, VE-20 LAB3, LABPETTE DIGITAL, VE-20, ch: SC, vol:0.2-20						
LABNET, VE-200 LAB3, LABPETTE DIGITAL, VE-200, ch: SC, vol:20-200					yes	
MATRIX, , MULTI-ELECTRAPETTE, , ch: MC, E, vol:0-25					yes	



METTLER-TOLEDO , VOLUMATE , , ch: SC, vol:0.1-2.5	yes					
METTLER-TOLEDO , VOLUMATE , , ch: SC, vol:0.5-10	yes	yes	yes	yes		
Pipetter Description	5984D	5984C	5984F	5984H	5984L	5984N
METTLER-TOLEDO , VOLUMATE , , ch: SC, vol:10-100						
METTLER-TOLEDO , VOLUMATE , , ch: SC, vol:20-200					yes	
METTLER-TOLEDO , VOLUMATE , , ch: SC, vol:2-20						
NICHIRYO , NICHIPET EX , , ch: SC, vol:10-100						
NICHIRYO , NICHIPET EX , , ch: SC, vol:20-200						
NICHIRYO , NICHIPET EX , , ch: SC, vol:2-20						
NICHIRYO, 5000 , , , ch: SC, vol:10-50						
NICHIRYO, 5000 , , , ch: SC, vol:200-1000						yes
NICHIRYO, 5000 , , , ch: SC, vol:40-200					yes	
NICHIRYO, 5000F, FIXED VOLUME , , ch: SC, vol:10-10	yes	yes	yes	yes		
NICHIRYO, 5000F, FIXED VOLUME , , ch: SC, vol:200-200						
NICHIRYO, 5000V-L , , , ch: SC, vol:40-200						
NICHIRYO, 5000V-M, TRIPLE VOLUME, ch: SC, vol:200-1000200/500/1000µl						yes
NICHIRYO, 5000V-M, TRIPLE VOLUME, ch: SC, vol:50-20050/100/200µl					yes	
NICHIRYO, 7000, NICHIPET DIGITAL , , ch: MC-8, vol:40-200					yes	
NICHIRYO, 7000, NICHIPET DIGITAL , , ch: MC-8, vol:40-200						
OXFORD®, , BENCHMATE®, , ch: MC-12, vol:40-200					yes	
OXFORD®, , BENCHMATE®, , ch: SC, vol:0.5-10	yes	yes	yes	yes		
OXFORD®, , BENCHMATE®, , ch: SC, vol:10-50						
OXFORD®, , BENCHMATE®, , ch: SC, vol:200-1000						yes
OXFORD®, , BENCHMATE®, , ch: SC, vol:40-200					yes	
RAININ , , PIPET-PLUS, R-10, ch: SC, vol:0.5-10	yes	yes				
RAININ , , PIPET-PLUS, R-100, ch: SC, vol:10-100					yes	
RAININ , , PIPET-PLUS, R-1000, ch: SC, vol:100-1000						yes



RAININ, , PIPET-PLUS, R-2, ch: SC, vol:0.1-2	yes	yes				
RAININ, , PIPET-PLUS, R-20, ch: SC, vol:2-20			yes	yes		
Pipetter Description	5984D	5984C	5984F	5984H	5984L	5984N
RAININ, , PIPET-PLUS, R-200, ch: SC, vol:50-200					yes	
SOCOREX, 811, , , ch: SC, vol:100-100					yes	
SOCOREX, 811, , , ch: SC, vol:5-5			yes			
SOCOREX, 815, Acura Fixed Volume, , ch: SC, vol:200-200						
SOCOREX, 821, ch: SC, vol:200-1000						yes
SOCOREX, 821, , , ch: SC, vol:50-200					yes	
SOCOREX, 821, , , ch: SC, vol:5-50						
SOCOREX, 822, CALIBRA®, , ch: SC, vol:0.4-2			yes			
SOCOREX, 822, CALIBRA®, , ch: SC, vol:100-1000						yes
SOCOREX, 822, CALIBRA®, , ch: SC, vol:10-100						
SOCOREX, 822, CALIBRA®, , ch: SC, vol:1-10	yes	yes	yes	yes		
SOCOREX, 822, CALIBRA®, , ch: SC, vol:1-10			yes			
SOCOREX, 822, CALIBRA®, , ch: SC, vol:20-200					yes	
SOCOREX, 822, CALIBRA®, , ch: SC, vol:2-20						
SOCOREX, 825, Acura, , ch: SC, vol:0.1-2	yes	yes	yes	yes		
SOCOREX, 825, Acura, , ch: SC, vol:0.5-10	yes	yes	yes	yes		
SOCOREX, 825, Acura, , ch: SC, vol:100-1000						yes
SOCOREX, 825, Acura, , ch: SC, vol:10-100					yes	
SOCOREX, 825, Acura, , ch: SC, vol:20-200					yes	
SOCOREX, 825, Acura, , ch: SC, vol:2-20						
SOCOREX, 851, ACURA ANALOGIC, ch: MC-8, vol:50-200					yes	
SOCOREX, 851, ACURA ANALOGIC, , ch: MC-8,12, vol:5-50						
SOCOREX, 852, CALIBRA®, , ch: MC-8, vol:10-100						
SOCOREX, 852, CALIBRA®/WHEATON, , ch: MC-8, vol:1-10	yes	yes				



TITERTEK®/ PLUS ICN , ch: MC-8, E, vol:5-100					yes	
TITERTEK® / PLUS, ch: MC-12, vol:5-50						
Pipetter Description	5984D	5984C	5984F	5984H	5984L	5984N
TITERTEK®, DIGITAL, ch: MC-12, vol:50-300					yes	
TITERTEK®, DIGITAL, ch: MC-12, vol:5-50						
TITERTEK®,DIGITAL, ch: SC, vol:200-1000						yes
TITERTEK®, DIGITAL, ch: SC, vol:40-200					yes	
TITERTEK®,DIGITAL, ch: SC, vol:5-40					yes	
VISTALAB,MLA®, ch: SC, vol:100-100, 1025					yes	
VISTALAB , MLA®, , ch: SC, vol:20-20					yes	
VISTALAB , MLA®, , ch: SC, vol:20-20						
VISTALAB , MLA®, , ch: SC, vol:50-200						
VISTALAB , MLA®, , ch: SC, vol:50-50					yes	
VISTALAB,MLA®, Precision D Calibratable, ch: SC, vol:10-10					yes	
VISTALAB,MLA®, Precision D Calibratable, ch:SC, vol:200-200						
VISTALAB,MLA®, Precision D Calibratable, ch:SC, vol:25-25					yes	
VISTALAB,MLA®, Precision D Calibratable, ch: SC, vol:5-5					yes	
VISTALAB,OVATION, ch: MC-12, vol:5-250						
VISTALAB, OVATION, ch: SC, vol:0.2-2	yes					
VISTALAB,OVATION,ch: SC, vol:10-100					yes	
VISTALAB,OVATION,ch: SC, vol:1-10		yes	yes	yes		
VISTALAB,OVATION,ch: SC, vol:2-20		yes	yes	yes		
VISTALAB , OVATION, gray cone installed, ch: SC, vol:20-200						
VISTALAB , , OVATION, OVATION, ch: MC-8, vol:5-250						
VISTALAB , , OVATION, white cone installed, ch: SC, vol:100-1000						yes
VISTALAB , , OVATION-white cone, old version, white cone, ch: SC, vol:20-200						
VWR Signature, , International, , ch: MC 12 E, vol:10-200						



VWR Signature, , International, , ch: MC 12, vol:50-300						
VWR Signature, , International, , ch: MC 12, vol:5-50						
Pipetter Description	5984D	5984C	5984F	5984H	5984L	5984N
VWR Signature, , International, , ch: MC 8, vol:1-10	yes	yes	yes	yes		
VWR Signature, , International, , ch: SC, vol:0.1-2	yes	yes	yes	yes		
VWR Signature, , International, , ch: SC, vol:0.5-10	yes	yes				
VWR Signature, , International, , ch: SC, vol:10-100					yes	
VWR Signature, , International, , ch: SC, vol:20-200					yes	
VWR Signature, , International, , ch: SC, vol:2-20						